

Custom Solution - 6 components in HNO3 5%

Art. ID LM24-CUS-72552
Unit 2 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	1000	µg/mL		
	Magnesium (Mg)	[7439-95-4]	2000	µg/mL		
	Potassium (K)	[7440-09-7]	10000	µg/mL		
	Sodium (Na)	[7440-23-5]	3000	µg/mL		
	Lithium (Li)	[7439-93-2]	100	µg/mL		
	Phosphorus (P)	[7723-14-0]	1000	µg/mL		