

Custom Solution - 11 components 10 mg/L each of in HCl 10%

Art. ID LM24-CUS-72555
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Selenium (Se)	[7782-49-2]	10	mg/L		
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Tin (Sn)	[7440-31-5]	10	mg/L		
	Antimony (Sb)	[7440-36-0]	10	mg/L		
	Mercury (Hg)	[7439-97-6]	10	mg/L		
	Titanium (Ti)	[7440-32-6]	10	mg/L		
	Tungsten (W)	[7440-33-7]	10	mg/L		
	Iridium (Ir)	[7439-88-5]	10	mg/L		
	Gold (Au)	[7440-57-5]	10	mg/L		
	Ruthenium (Ru)	[7440-18-8]	10	mg/L		
	Platinum (Pt)	[7440-06-4]	10	mg/L		