

### **Custom Solution - 2 components 100 mg/L each of in Methanol**

Art. ID LM24-CUS-72559  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bisphenol A D16	[96210-87-6]	100	mg/L		
	4-n-Nonylphenol D8 [358 730-95-7]	[358730-95-7]	100	mg/L		