

Custom Solution - 39 components in Methanol + Extended Shelf-Life (Option D): + 12 months

Art. ID LM24-CUS-72577
Unit 5 x 1 mL ampoule
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months + Extended Shelf-Life (Option D) + 12 months. (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,1-Dichloroethene	[75-35-4]	2000	µg/mL		
	Methylene chloride (Dichloromethane)	[75-09-2]	20000	µg/mL		
	trans-1,2-Dichloroethene	[156-60-5]	2000	µg/mL		
	1,1-Dichloroethane	[75-34-3]	2000	µg/mL		
	cis-1,2-Dichloroethene	[156-59-2]	2000	µg/mL		
	Chloroform (Trichloromethane)	[67-66-3]	2000	µg/mL		
	1,1,1-Trichloroethane	[71-55-6]	2000	µg/mL		
	Trichloroethene	[79-01-6]	2000	µg/mL		
	Tetrachloroethene	[127-18-4]	2000	µg/mL		
	Carbon tetrachloride (Tetrachloromethane)	[56-23-5]	2000	µg/mL		
	1,2-Dichloroethane	[107-06-2]	2000	µg/mL		
	1,2-Dichloropropane	[78-87-5]	2000	µg/mL		
	Bromodichloromethane	[75-27-4]	2000	µg/mL		
	1,1,2-Trichloroethane	[79-00-5]	2000	µg/mL		
	1,1,2,2-Tetrachloroethane	[79-34-5]	2000	µg/mL		
	Dibromochloromethane	[124-48-1]	2000	µg/mL		
	1,1,1,2-Tetrachloroethane	[630-20-6]	2000	µg/mL		
	1,2-Dibromoethane	[106-93-4]	2000	µg/mL		
	Bromoform (Tribromomethane)	[75-25-2]	2000	µg/mL		
	1,2,3-Trichloropropane	[96-18-4]	2000	µg/mL		
	Hexachloro-1,3-butadiene	[87-68-3]	2000	µg/mL		

Toluene	[108-88-3]	2000	µg/mL
Styrene	[100-42-5]	2000	µg/mL
p-Xylene	[106-42-3]	2000	µg/mL
m-Xylene	[108-38-3]	2000	µg/mL
o-Xylene	[95-47-6]	2000	µg/mL
Naphthalene	[91-20-3]	2000	µg/mL
Isopropylbenzene (Cumol)	[98-82-8]	2000	µg/mL
Ethylbenzene	[100-41-4]	2000	µg/mL
Chlorobenzene	[108-90-7]	2000	µg/mL
Benzene	[71-43-2]	2000	µg/mL
1,4-Dichlorobenzene	[106-46-7]	2000	µg/mL
1,3-Dichlorobenzene	[541-73-1]	2000	µg/mL
1,2-Dichlorobenzene	[95-50-1]	2000	µg/mL
1,2,4-Trichlorobenzene	[120-82-1]	2000	µg/mL
Methyl-tert.butylether (MTBE)	[1634-04-4]	2000	µg/mL
Ethyl-tert-butylether (ETBE)	[637-92-3]	2000	µg/mL
1,2,3-Trichlorobenzene	[87-61-6]	2000	µg/mL
1,3,5-Trichlorobenzene	[108-70-3]	2000	µg/mL