

### **Custom Solution - Kit of 9 components 100 mg/L each of in Acetonitrile**

Art. ID LM24-CUS-72603  
Unit 3 x 1 mL ampoule  
Deliverydetails De Minimis UN1648 II E2

#### Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Trinexapac (free acid)	[143294-89-7]	100	mg/L		
	Trinexapac-ethyl	[95266-40-3]	100	mg/L		
	Triclopyr	[55335-06-3]	100	mg/L		
	Metobromuron	[3060-89-7]	100	mg/L		
	Dichlorprop-P	[15165-67-0]	100	mg/L		
	Diclofop-methyl	[51338-27-3]	100	mg/L		
	Cyhalofop-butyl	[122008-85-9]	100	mg/L		
	Thiencarbazone-methyl	[317815-83-1]	100	mg/L		
	2-Hydroxy-propoxycarbaz one	[496925-02-1]	100	mg/L		