

Custom Solution - 6 components in HNO3 5%

Art. ID LM24-CUS-72604
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	250	mg/L		
	Antimony (Sb)	[7440-36-0]	250	mg/L		
	Silicon (Si)	[7440-21-3]	250	mg/L		
	Tin (Sn)	[7440-31-5]	50	mg/L		
	Tellurium (Te)	[13494-80-9]	250	mg/L		
	Titanium (Ti)	[7440-32-6]	25	mg/L		