

Custom Solution - 26 components in HNO3 5%

Art. ID	LM24-CUS-72605
Unit	1 x 250 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	250	mg/L		
	Arsenic (As)	[7440-38-2]	75	mg/L		
	Boron (B)	[7440-42-8]	75	mg/L		
	Barium (Ba)	[7440-39-3]	75	mg/L		
	Beryllium (Be)	[7440-41-7]	2,5	mg/L		
	Bismuth (Bi)	[7440-69-9]	250	mg/L		
	Calcium (Ca)	[7440-70-2]	2500	mg/L		
	Cadmium (Cd)	[7440-43-9]	25	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	50	mg/L		
	Copper (Cu)	[7440-50-8]	250	mg/L		
	Iron (Fe)	[7439-89-6]	500	mg/L		
	Mercury (Hg)	[7439-97-6]	50	mg/L		
	Potassium (K)	[7440-09-7]	2500	mg/L		
	Lithium (Li)	[7439-93-2]	75	mg/L		
	Magnesium (Mg)	[7439-95-4]	2500	mg/L		
	Manganese (Mn)	[7439-96-5]	125	mg/L		
	Sodium (Na)	[7440-23-5]	2500	mg/L		
	Nickel (Ni)	[7440-02-0]	125	mg/L		
	Phosphorus (P)	[7723-14-0]	250	mg/L		
	Lead (Pb)	[7439-92-1]	125	mg/L		
	Selenium (Se)	[7782-49-2]	250	mg/L		
	Strontium (Sr)	[7440-24-6]	125	mg/L		
	Thallium (Tl)	[7440-28-0]	125	mg/L		
	Vanadium (V)	[7440-62-2]	125	mg/L		
	Zinc (Zn)	[7440-66-6]	250	mg/L		