

1,3-Bis(isopropylamino)propan-2-ol Dihydrochloride (1,3-bis(isopropylamino)propan-2-ol Dihydrochloride)

Art. ID USP-1A05580

Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,3-Bis[(1-methylethyl) amino]-2-propanol Dihydrochloride	[73313-36-7]				