

Custom Solution - 23 components in HNO3 5%

Art. ID	LM24-CUS-72608
Unit	1 x 100 mL bottle
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium (Li)	[7439-93-2]	5	mg/L		
	Beryllium (Be)	[7440-41-7]	5	mg/L		
	Vanadium (V)	[7440-62-2]	5	mg/L		
	Chromium (Cr)	[7440-47-3]	5	mg/L		
	Cobalt (Co)	[7440-48-4]	5	mg/L		
	Nickel (Ni)	[7440-02-0]	5	mg/L		
	Molybdenum (Mo)	[7439-98-7]	5	mg/L		
	Silver (Ag)	[7440-22-4]	5	mg/L		
	Cadmium (Cd)	[7440-43-9]	5	mg/L		
	Tin (Sn)	[7440-31-5]	5	mg/L		
	Thallium (Tl)	[7440-28-0]	5	mg/L		
	Lead (Pb)	[7439-92-1]	5	mg/L		
	Uranium (U)	[7440-61-1]	5	mg/L		
	Barium (Ba)	[7440-39-3]	20	mg/L		
	Boron (B)	[7440-42-8]	125	mg/L		
	Titanium (Ti)	[7440-32-6]	125	mg/L		
	Aluminium (Al)	[7429-90-5]	125	mg/L		
	Manganese (Mn)	[7439-96-5]	125	mg/L		
	Iron (Fe)	[7439-89-6]	125	mg/L		
	Copper (Cu)	[7440-50-8]	125	mg/L		
	Zinc (Zn)	[7440-66-6]	125	mg/L		
	Strontium (Sr)	[7440-24-6]	125	mg/L		
	Tellurium (Te)	[13494-80-9]	125	mg/L		