

Custom Solution - 4 components in HNO3 5%

Art. ID LM24-CUS-72611
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	5	mg/L		
	Selenium (Se)	[7782-49-2]	5	mg/L		
	Antimony (Sb)	[7440-36-0]	5	mg/L		
	Mercury (Hg)	[7439-97-6]	1	mg/L		