

Custom Solution - 2 components 1000 mg/L each of in HNO3 5%, HF traces

Art. ID LM24-CUS-72628
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Platinum (Pt)	[7440-06-4]	1000	mg/L		
	Tellurium (Te)	[13494-80-9]	1000	mg/L		