

Custom Solution - Kit of 8 components 0.2 mg/L each of in Acetonitrile

Art. ID	LM24-CUS-72631
Unit	5 x 1 mL ampoule
Deliverydetails	De Minimis UN1648 II E2

Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Monuron	[150-68-5]	0,2	mg/L		
	Dimethachlor metabolite CGA 369873 sodium	[2387071-47-6]	0,2	mg/L		
5-(Trifluoromethyl)-2-pyridone	2-Hydroxy-5-(trifluoromethyl)pyridine	[33252-63-0]	0,2	mg/L		
	Metaldehyde	[108-62-3]	0,2	mg/L		
	Imazalil	[35554-44-0]	0,2	mg/L		
	Terbutylazine Metabolite CGA 324007	[309923-18-0]	0,2	mg/L		
	6-(Tert-Butylamino)-3-Methyl-1,3,5-Triazine-2,4(1H,3H)-Dione	[2206682-85-9]	0,2	mg/L		
	Dimethachlor Metabolite CGA 373464	[1196157-87-5]	0,2	mg/L		