

### **Custom Solution - Terbutryn [886-50-0] 1000 mg/L in Methanol**

Art. ID LM24-CUS-72634  
Unit 1 x 1 ampule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation 12 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
|                  | Terbutryn         | [886-50-0] | 1000                | mg/L |        |        |