

Custom Solution - 17 components in HNO3 2%

Art. ID LM24-CUS-72643
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10	mg/L		
	Boron (B)	[7440-42-8]	75	mg/L		
	Barium (Ba)	[7440-39-3]	1	mg/L		
	Calcium (Ca)	[7440-70-2]	50	mg/L		
	Cadmium (Cd)	[7440-43-9]	1	mg/L		
	Cobalt (Co)	[7440-48-4]	5	mg/L		
	Chromium (Cr)	[7440-47-3]	1	mg/L		
	Copper (Cu)	[7440-50-8]	5	mg/L		
	Iron (Fe)	[7439-89-6]	5	mg/L		
	Magnesium (Mg)	[7439-95-4]	6	mg/L		
	Manganese (Mn)	[7439-96-5]	1	mg/L		
	Nickel (Ni)	[7440-02-0]	5	mg/L		
	Lead (Pb)	[7439-92-1]	5	mg/L		
	Tin (Sn)	[7440-31-5]	2	mg/L		
	Titanium (Ti)	[7440-32-6]	5	mg/L		
	Vanadium (V)	[7440-62-2]	2	mg/L		
	Zinc (Zn)	[7440-66-6]	5	mg/L		