

Custom Solution - 8 components in HNO3 5%

Art. ID LM24-CUS-71250-2X100ML
Unit 2 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 12 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1000	µg/mL		
	Calcium (Ca)	[7440-70-2]	1000	µg/mL		
	Iron (Fe)	[7439-89-6]	1000	µg/mL		
	Magnesium (Mg)	[7439-95-4]	1000	µg/mL		
	Sodium (Na)	[7440-23-5]	1000	µg/mL		
	Potassium (K)	[7440-09-7]	1000	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	50	µg/mL		
	Titanium (Ti)	[7440-32-6]	50	µg/mL		