

Custom Solution - 4 components in HNO3 0.2%

Art. ID LM24-CUS-72648
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium (Na+)		1000	mg/L		
	Potassium (K+)		200	mg/L		
	Magnesium (Mg ²⁺)		2000	mg/L		
	Calcium (Ca ²⁺)		5000	mg/L		