

**Custom Solution - 5 components 100 mg/L each of in HNO3 5%**

Art. ID LM24-CUS-71789-2X100ML  
Unit 2 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

**Description**

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Strontium (Sr)	[7440-24-6]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		