

Custom solution - 5 components 10 mg/L each of in 2% HNO₃

Art. ID LM24-CUS-72659
Unit 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	10	mg/L		
	Bismuth (Bi)	[7440-69-9]	10	mg/L		
	Indium (In)	[7440-74-6]	10	mg/L		
	Iridium (Ir)	[7439-88-5]	10	mg/L		
	Rhodium (Rh)	[7440-16-6]	10	mg/L		