

Custom Solution - 9 components 100 mg/L each of in HCl 15%

Art. ID LM24-CUS-72669
Unit 1 x 100 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1789 II E2

Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Germanium (Ge)	[7440-56-4]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		
	Sulfur (S)	[7704-34-9]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Tungsten (W)	[7440-33-7]	100	mg/L		
	Zirconium (Zr)	[7440-67-7]	100	mg/L		