

Cast Iron, CRM, Ductile (Nodular) Iron, Block 48 mm x 42 mm x 12 mm

Art. ID BAS-SCRM-670/25
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,178 ± 0,024	%		
	Cerium (Ce)	[7440-45-1]	0,0092 ± 0,0014	%		
	Chromium (Cr)	[7440-47-3]	0,488 ± 0,011	%		
	Copper (Cu)	[7440-50-8]	0,978 ± 0,019	%		
	Magnesium (Mg)	[7439-95-4]	0,0414 ± 0,0021	%		
	Manganese (Mn)	[7439-96-5]	0,324 ± 0,008	%		
	Molybdenum (Mo)	[7439-98-7]	0,0110 ± 0,0014	%		
	Nickel (Ni)	[7440-02-0]	0,965 ± 0,010	%		
	Silicon (Si)	[7440-21-3]	2,152 ± 0,017	%		
	Titanium (Ti)	[7440-32-6]	0,0820 ± 0,0030	%		
	Vanadium (V)	[7440-62-2]	0,0125 ± 0,0006	%		
	Sulfur (S)	[7704-34-9]	0,0085 ± 0,0012	%		
	Phosphorus (P)	[7723-14-0]	0,0466 ± 0,0020	%		