

rac-Itraconazole Isopropyl Analog

(4-(4-{4-[4-((2RS,4SR)-2-[(1H-1,2,4-Triazol-1-yl)methyl]-2-(2,4-dichlorophenyl)-1,3-dioxolan-4-yl)methoxy]phenyl}piperazin-1-yl)phenyl)-2-isopropyl-2H-1,2,4-triazol-3(4H)-one)

Art. ID USP-1A08040

Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	rac-Itraconazole Isopro pyl Analog	[89848-49-7]				