

**Custom Solution - 1 component, Cu 100 µg/L in HNO<sub>3</sub> 2%**

Art. ID LM24-CUS-72685  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

**Description**

Conservation 8 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	100	µg/L		