

Cinchonidine Sulfate Dihydrate
(R-Quinolin-4-yl[(2S,4S,5R)-5-vinylquinuclidin-2-yl]methanol sulfate dihydrate)

Art. ID USP-1A08690
Unit 25 mg

Description

Note: PAI products are different from official USP reference standards. These products are released using a separate process developed by USP's subject matter experts. The release process follows rigorous internal policies, standard operating procedures, and USP Quality Management System requirements. USP is an ISO9001:2015 certified facility.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	R-Quinolin-4-yl[(2S,4S,5R)-5-vinylquinuclidin-2-yl]methanol sulfate dihydrate	[n/a]				