

### **Custom Solution - Acrylamide-d3 [122775-19-3] 1000 mg/L in Methanol**

Art. ID LM24-CUS-70794-2x1ML  
Unit 2 x 1 mL ampoule  
Deliverydetails De Minimis UN1230 (6.1) II E2

#### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acrylamide-d3	[122775-19-3]	1000	mg/L		