

25-Dihydroxyvitamin D3 (25,26,27-13C3, 98%) 50 µg/mL in Ethanol CP 95%

Art. ID CIL-CLM-11583-B
Unit 1 mL
Deliverydetails Dry ice shipment UN1170 III E1

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|----------------------------|---------------------|-------|--------|--------|
| | 25-Dihydroxyvitamin D3 (25,26,27-13C3, 98%) | [19356-17-3]/un labeled | 50 | µg/mL | | |