

**3'-Hydroxycotinine (mix of stereoisomers) (2',3',4'-13C3, 99%) 100 µg/mL in Methanol CP 95%**

Art. ID CIL-CLM-11520-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3'-Hydroxycotinine (mix of stereoisomers) (2',3',4'-13C3, 99%)	[79581-34-3]/un labeled	100	µg/mL		