

Bisphenol AF (ring-13C12, 99%) 1 mg/mL in Methanol

Art. ID CIL-CLM-9776-10X-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bisphenol AF (ring-13C1 2,99%)	[1478-61-1]/unl abelled	1	mg/mL		