

Bis(2-ethylhexyl)phosphate (D34, 97%) 100 µg/mL in Acetonitrile CP 97%

Art. ID CIL-DLM-11382-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------------|----------------|---------------------|-------|--------|--------|
| | Bis(2-ethylhexyl) Phosp hate-d34 | [1773493-20-1] | 100 | µg/mL | | |