

**Bis(2-ethylhexyl)phosphate (D34, 97%) 100 µg/mL in Acetonitrile CP 97%**

Art. ID CIL-DLM-11382-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bis(2-ethylhexyl) Phosp hate-d34	[1773493-20-1]	100	µg/mL		