

Bisphenol A (ring-3,3',5,5'-D4, 97%) 1000 µg/mL in Methanol-OD

Art. ID CIL-DLM-2774-10X-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bisphenol A-3,3',5,5'-D 4	[347841-41-2]	1000	µg/mL		