

### **Bisphenol AF (unlabeled) 1000 µg/mL in Methanol**

Art. ID CIL-ULM-9779-10X-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

#### Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bisphenol AF	[1478-61-1]	1000	µg/mL		