

## **Custom Solution - Phosphates (PO43-) 10000 mg/L in Water**

Art. ID LM24-CUS-70660-3X100ML  
Unit 3 x 100 mL  
Deliverydetails No Dangerous Good /not restricted

### Description

Conservation: 24 months (originally closed bottle).

| Text/Information | Analyte/Parameter  | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------|------------|---------------------|------|--------|--------|
|                  | Phosphates (PO43-) |            | 10000               | mg/L |        |        |