

### **Custom Solution - 15 components 100 mg/L each of in HNO3 2%**

Art. ID LM24-CUS-72747-LOT A  
Unit 15 x 500 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### Description

Conservation 24 months (in originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Lithium (Li)	[7439-93-2]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		
	Sulfur (S)	[7704-34-9]	100	mg/L		