

## **Custom Solution - 8 components 2000 mg/L each of in Ethyl Acetate**

Art. ID	LM24-CUS-72751
Unit	3 x 1 mL ampoule
Deliverydetails	De Minimis UN1173 II E2

### Description

Conservation 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phthalic acid, bis-ethyl ester (Bis-ethyl phthalate, Diethyl phthalate)	[84-66-2]	2000	mg/L		
	Phthalic acid, bis-butyl ester (Bis-butyl phthalate, Di-n-butyl phthalate)	[84-74-2]	2000	mg/L		
	Phthalic acid, bis-isobutyl ester (Bis-isobutyl phthalate)	[84-69-5]	2000	mg/L		
	Phthalic acid, benzylbutyl ester (Benzylbutyl phthalate)	[85-68-7]	2000	mg/L		
	Adipic acid-bis-2-ethylhexyl ester (Bis(2-ethylhexyl)adipate)	[103-23-1]	2000	mg/L		
	Phthalic acid, bis-2-ethylhexyl ester (Bis-2-ethylhexyl phthalate, Bis(2-ethylhexyl)phthalate)	[117-81-7]	2000	mg/L		
	Phthalic acid, bis-1-octyl ester (Bis-1-octyl phthalate, Di-n-octylphthalate)	[117-84-0]	2000	mg/L		
	Phthalic acid, bis-isononyl ester (Diisononyl Phthalate) (mixture of	[28553-12-0]	2000	mg/L		

branched chain isomers

)