

## **Custom Solution - 5 components in HNO3 5% for ICP-MS**

Art. ID LM24-CUS-72761  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

### Description

Conservation 24 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium (Na)	[7440-23-5]	2500	mg/L		
	Magnesium (Mg)	[7439-95-4]	2500	mg/L		
	Potassium (K)	[7440-09-7]	500	mg/L		
	Calcium (Ca)	[7440-70-2]	2500	mg/L		
	Silicon (Si)	[7440-21-3]	500	mg/L		