

## **2,6-Dinitrotoluene-d3 100 µg/mL in Methanol**

Art. ID LM24-N-13400-4956D-100ME1  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,6-Dinitrotoluene-alph a,alpha,alpha-d3	[93951-90-7]	100	µg/mL		