

1,3,5-Trimethylbenzene-2,4,6-d3 100 µg/mL in Methanol

Art. ID LM24-N-13400-5008-100ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	D3-1,3,5-Trimethylbenzene	[38574-14-0]	100	µg/mL		