

## **Propionaldehyde-DNPH-d3 10 µg/mL in Acetonitrile**

Art. ID LM24-N-13430-0217D-10AN5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Propionaldehyde 2,4-Dinitrophenylhydrazone-d3	[259824-53-8]	10	µg/mL		