

Propionaldehyde-DNPH-d3 10 µg/mL in Acetonitrile (based on active substance)

Art. ID LM24-N-13430-0217D-W10AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
based on active substan ce	Propionaldehyde 2,4-Din itrophenylhydrazone-d3	[259824-53-8]	10	µg/mL		