

(-)-2-Methylisoborneol-d3 100 µg/mL in Methanol

Art. ID LM24-N-13500-0632D-100ME1
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------------|---------------|---------------------|-------|--------|--------|
| | (-)-2-Methyl Isoborneol -d3 | [135441-89-3] | 100 | µg/mL | | |