

Tetrafluoroterephthalic acid 1000 µg/mL in Methanol

Art. ID LM24-N-13700-0335-1000ME5

Unit 5 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tetrafluoroterephthalic acid	[652-36-8]	1000	µg/mL		