

## **Sodium Dodecyl-sulfate-d25 100 µg/mL in Methanol**

Art. ID LM24-N-13700-0338D-100ME5  
Unit 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

### Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sodium dodecyl sulfate- d25	[110863-24-6]	100	µg/mL		