

Trifluoroacetic acid-d1 1000 µg/mL in Acetonitrile

Art. ID LM24-N-13700-0344D-1000AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Trifluoroacetic acid-d	[599-00-8]	1000	µg/mL		