

Captan-3,3,4,5,6,6-d6 10 µg/mL in Acetonitrile

Art. ID LM24-N-15900-0611D-10AN5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Captan-3,3,4,5,6,6-D6	[1330190-00-5]	10	µg/mL		