

### **Sulcotrione-d7 100 µg/mL in Acetonitrile**

Art. ID LM24-N-15900-5742D-100AN1  
Unit 1 mL  
Deliverydetails De Minimis UN1648 II E2

#### Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|------------|---------------------|-------|--------|--------|
|                  | Sulcotrione-d7    |            | 100                 | µg/mL |        |        |