

Alachlor-ESA sodium salt-d3 100 µg/mL in Methanol

Art. ID LM24-N-15900-5794-100ME1
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------|------------|---------------------|-------|--------|--------|
| | Alachlor-ESA sodium sal t-d3 | [n/a] | 100 | µg/mL | | |