

**Alachlor-ESA sodium salt-d3 100 µg/mL in Methanol (based on active substance)**

Art. ID LM24-N-15900-5794-W100ME1  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
based on active substance	Alachlor-ESA sodium salt-d3	[n/a]	100	µg/mL		