

Alachlor-ESA sodium salt-d3 10 µg/mL in Methanol (based on active substance)

Art. ID LM24-N-15900-5794-W10ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|---------------------------|-----------------------------|------------|---------------------|-------|--------|--------|
| based on active substance | Alachlor-ESA sodium salt-d3 | [n/a] | 10 | µg/mL | | |