

Alachlor-ESA sodium salt-d3 10 µg/mL in Methanol (based on active substance)

Art. ID LM24-N-15900-5794-W10ME5
Unit 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
based on active substance	Alachlor-ESA sodium salt-d3	[n/a]	10	µg/mL		