

2,4-D-sodium-monohydrate 10 µg/mL in Acetonitrile

Art. ID LM24-N-15900-5806-10AN10
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4-D-sodium salt monohydrate	[2702-72-9]	10	µg/mL		